Brainy

SMART THERMOMETER

INFRARED
AGU SHE7

- DE BEDIENUNGSANLEITUNG
Fassung 2.0
- FR MODE D'EMPLOI
Version 2.0
- EN USER MANUAL
Version 2.0
- RU РУКОВОДСТВО
ПО ЭКСПЛУАТАЦИИ
Версия 2.0

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CONTENTS

1	Introduction	51
2	List of symbols	53
3	Scope of application	54
4	Complete set	54
5	Specifications	54
6	Device description	56
7	Safety instructions	58
8	Preparation for use	60
9	Operating principles and procedures	62
10	Device cleaning	66
11	Troubleshooting	67
12	Storage, transportation and operation terms	68
13	Disposal	68
14	Certification	69
15	Manufacturer's warranty	71

1 INTRODUCTION

Dear friends, thank you for choosing our products!

Smart infrared thermometer **AGU SHE7** is an advanced high-quality product that is compliant with international standards.

Smart infrared thermometer **AGU SHE7** has been designed for body temperature measurement.

It performs high-precision temperature measurements through the use of unique technology. Power-on self-test of the system ensures the required measuring accuracy.

This device is suitable for people of all ages.

Benefits of smart infrared thermometer **AGU SHE7**:

Multiple purpose (wide measurement range).

The device has a wide measurement range from 0°C to 100°C, which allows it to be used as a forehead or ear thermometer, as well as for measuring:

- liquid surface temperature in a baby bottle (for instance, milk);
- water surface temperature in a baby bath;
- ambient temperature and solid object surface temperature.

Measuring temperature in a just few seconds.

Innovative technology with the use of infrared sensor allows measuring temperature in the external ear canal in just a 1 second.

Accuracy and reliability.

The unique design of the device with built-in infrared sensor provides accurate and reliable measurement results.

Ease of use.

Smart infrared thermometer **AGU SHE7** allows conducting temperature measurements easily and simply.

The temperature can be measured even in a sleeping child without disturbing.

The measurement is fast, which is very convenient when taking a child's temperature.

Automatic display of readings in the device memory.

The device automatically displays the last measurement for 2 seconds after switching on.

Reproduction of the last few results.

You can recall the last 10 measurement results, which makes it possible to more effectively monitor temperature changes.

Data storage and viewing in the application due to wireless transmission.

Having installed the AGU application on your smartphone, you can save and view temperature measurement statistics.

Please read this manual carefully before using the device.

2 LIST OF SYMBOLS

Symbol	Meaning
	The CE marking with the Registration Number of the Notified Body. This denotes the compliance of European Medical Device Directive 93/42/EEC.
	The product complies with the basic requirements of TR CU 020/2011 «Electromagnetic compatibility of technical means».
	WEEE (Waste Electrical and Electronic Equipment Directive). The symbol on the product or its package means that this product does not fall under the category of domestic waste. To avoid possible damage to the environment and human health, separate such wastes from others and dispose of them in accordance with accepted standards.
	Device classification type BF.
IP22	Ingress protection rating. Leading digit (protection against ingress of solid foreign objects): 2 – protection against ingress of solid objects more than 12 mm in size; fingers or other objects with a maximum length of 80 mm, or solid objects. Second digit (protection against ingress of foreign liquids): 2 – protection against vertically falling water drops and objects when enclosure tilted up to 15° (normally positioned equipment).
	Warning/Attention.
	Read the instruction manual before use.

3 SCOPE OF APPLICATION

Smart infrared thermometer **AGU SHE7** is intended for measuring body temperature on the forehead and in the external ear canal in the home environment. Do not use the device for purposes other than intended. This thermometer can be used for children and adults.

4 COMPLETE SET


1. Smart non-contact digital infrared thermometer **AGU SHE7** – 1 pc.
2. Batteries 1.5 V AAA – 2 pcs.
3. Storage case – 1 pc.
4. Instruction manual – 1 pc.

5 SPECIFICATIONS


Type	Smart infrared thermometer
Model	AGU SHE7
Temperature measurement range	Human body: 34°C – 43°C (93.2°F – 109.4°F). Objects: 0°C – 100°C (32°F – 212°F)
Measuring accuracy	Human body: 35°C – 42°C (95°F – 107.6°F) ±0.2°C (0.4°F). Objects: ±2°C (±4°F) or ±5%
Operating environment temperature	15°C – 35°C (59°F – 95°F). Relative humidity: up to 95% (non-condensing)
Storage conditions	-25°C...+55°C (-13°F...+131°F). Relative humidity: up to 95% (non-condensing)
Transport conditions	-25°C...+55°C (-13°F...+131°F). Relative humidity: up to 85% (non-condensing)
Operating conditions	15°C – 35°C (59°F – 95°F). Relative humidity: up to 95% (non-condensing)
Display resolution	0.1°C or 0.1°F
Power source	2 batteries 1.5 V AAA

Memory	10 records
Weight (without batteries), g	69
Dimensions (LxWxH), mm	134.3×42.4×47.6
Housing material	ABS plastic
Measurement time	1 second (in the external ear canal)
Display	segment
Auto power-off	in 30 seconds

6 DEVICE DESCRIPTION

DESCRIPTION OF THE DISPLAY ELEMENTS

7 SAFETY INSTRUCTIONS

The device must be used for its intended purpose, in accordance with this manual. The manufacturer is not responsible for any damages caused by improper use.

- Any modification of the device is prohibited.
- The device is not waterproof! Under no circumstances should you immerse it in water or other liquids. When cleaning, follow the instructions contained in «Device cleaning» section.
- After cleaning the measuring sensor, take a five-minute break before the next measurement so that the sensor gains the initial operating temperature.
- Do not use the device if you think it is damaged or notice anything unusual in its operation.
- Never disassemble the thermometer.
- It is important to clean the ear canal before measurement as the ear wax can affect the measurement results (usually too low).
- During early stages of illness with a rapid rise in temperature to high values, the effect of «white hyperthermia» can be observed – a condition which can be characterized by the peripheral vessel constriction, when the skin becomes pale and remains cold. In such cases, temperature measurements on the forehead should not be carried out, because the skin temperature will be low. THEREFORE, YOU SHOULD SELECT «EAR» MODE TO MEASURE THE TEMPERATURE IN THE EXTERNAL EAR CANAL.
- If the measurement results are contrary to the patient's general state or suspiciously low, repeat the measurements every 15 minutes or measure the body temperature with another thermometer in the mouth or in the rectum.
- The device includes parts that require careful handling.
- Please refer to «Basic specifications» section for storage and operation conditions.
- Ensure that children do not use the device unsupervised: some parts are small enough to be swallowed. Do not allow children to play with the thermometer.
- Remove the batteries from the device if it will not be used for extended period of time.
- Do not use the device near strong electromagnetic field sources, including cellular phones and radio installations. Maintain the required distance from these sources during operation of the device.
- Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated. In case of contamination,

clean the measuring sensor lens, wait for a few minutes before taking the next measurement so that the sensor gains the initial operating temperature.

- Do not pocket the thermometer before use as it can be heated.
- Do not measure temperature during or straight after nursing.
- When measuring temperature in babies, it is ideal to place the baby on a flat surface and turn the head on one side so that the ear is directed upward.
- Keep your thermometer out of the reach of children! The device should not be used by physically, nervously or mentally disordered people (including children), or by the people who have insufficient experience and knowledge, unless they are supervised or instructed by the persons responsible for their safety. You should supervise the children and not allow them to play with the thermometer.
- The user should not take food, drink or do physical exercises immediately prior to or during temperature measurement.
- Before starting measurement, please remove any dirt, hair or sweat from the measuring area.
- The environmental temperature in the apartment may differ in various rooms. Before measuring temperature, the user and the thermometer should stay in the same room (premises) with constant environmental conditions for at least 30 minutes.
- Always measure the temperature in the same area, as the readings may vary in different areas.
- Wait for at least 30 seconds after 3-5 consecutive measurements to ensure accuracy of readings.


Protect the device from:

- extreme temperatures;
- shocks;
- contamination and dust;
- direct sunlight.

This device is not intended to be a substitute for professional medical advice.

8 PREPARATION FOR USE

Remove the device from the package. Check the battery charge level.

When the batteries get weak, the battery warning symbol appears.

It is still possible to make measurements, but the batteries should be replaced as soon as possible. If the batteries are completely discharged, the display shows «Lo» and a low battery indicator.


In this case, the batteries must be replaced.

⚠ NOTE

- It is recommended to remove the batteries if the thermometer will not be used for extended period of time.
- Keep the batteries out of the reach of children and do not expose them to heat.
- Batteries must be disposed of in accordance with local environmental and organization policies.
- Do not use rechargeable batteries.

Battery replacement

1. Loosen the battery lid screw with a Phillips screwdriver. Remove the lid.
2. Insert 2 new batteries 1.5 V AAA into the battery compartment, observing polarity. Always use only new batteries.
3. Reinstall the battery lid and tighten the screw.

Dispose of the used batteries in accordance with the applicable legal regulations. It is prohibited to dispose of the batteries in the normal household waste.

DEVICE SETTINGS

Date and time settings

Before the first use of the device and after each replacement of batteries, date and time are displayed in the format 1-1 00:00. Set the following parameters in turn: temperature measurement unit, date and time.

Make settings in steps:




measurement unit – year – month – day – hours – minutes

1. Press «EAR» button.
2. Press and hold «SET» button for 5 seconds to activate the setting mode.
3. «C» or «F» measurement unit will be displayed.
 - Press «SET» button to change the parameter.
 - Press «MEMORY» button to confirm the entry.
4. The year indicator will be displayed, set the current date using «SET» and «MEM» buttons.
5. Set in turn the month, day, hours and minutes using «SET» and «MEM» buttons.
6. After setting all the parameters, the «OFF» indicator will be displayed and the device will switch off automatically.

Reproduction of the previous readings

You can recall up to 10 measurement values stored in the device memory.

1. Press «EAR» button to switch on the thermometer.
2. Press «MEMORY» button. The last measurement value stored in memory is displayed first.
3. Press «MEMORY» button one more time to view the previously stored value.
 - When recording new measurement values (after the 10th value), old values are deleted automatically.

Enable/disable sound

Press the «EAR» button to turn on the device. Then press the «SET» button to enable or disable the sound alert.

9 OPERATING PRINCIPLES AND PROCEDURES

There are 2 options for the device operation:

Option 1: **without mobile application.**

Option 2: **via mobile application.**


OPERATION WITHOUT MOBILE APPLICATION

Smart infrared thermometer **AGU SHE7** is intended for measuring temperature on the forehead and in external ear canal.


1. MEASURING BODY TEMPERATURE ON THE FOREHEAD

The thermometer takes temperature by measuring infrared heat generated in the eyebrow area and surrounding tissues.

1. Remove any dirt, hair or sweat before placing the temperature sensor in the measurement area.
2. Press «EAR» button to switch on the device.
3. Put the sensor tip with the cap to the temple and, holding «HEAD» button, swipe it smoothly across the forehead to the other temple.

4. Release head button «HEAD». You will hear a beep, and the result of the body temperature will appear on the display.

2. MEASURING BODY TEMPERATURE IN THE EAR



ATTENTION

- Some people produce different readings in their left and right ear. To avoid temperature variations, always perform measurements in the same ear.
- The measurement must not be taken in an ear affected by inflammatory diseases, after possible ear injuries or in the healing period. In all of these cases, please consult your doctor.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- Accumulation of ear wax on the sensor can affect the measuring accuracy and result in the spread of infection among the thermometer users. It is extremely important for each measurement to observe the rules of hygiene and keep the sensor clean. To clean the sensor, follow the instructions contained in the «Device cleaning» section. Every time before starting measurement, please make sure that the measuring sensor lens is not contaminated.
- After cleaning the measuring sensor lens with a disinfectant solution, wait for a few minutes before taking the next measurement so that the sensor gains the initial operating temperature.

The thermometer takes temperature by measuring infrared heat generated by ear drum and surrounding tissues. Energy is captured with a lens and converted to a temperature value. Temperature measurement should be carried out by placing the sensor tip in the external ear canal above the ear drum, and not over the surface of nearby anatomical areas. Otherwise, the measurement values may be too high.

Body temperature measurement in the external ear canal ensures accurate readings due to the same temperature of blood flowing towards the ear drum and the brain's temperature regulation center – hypothalamus. Therefore, changes in body temperature are faster and more accurately recorded in the external ear canal than in any other body area.

1. Press «EAR» button to switch on the device.
2. Remove the cap.

3. Make sure the ear canal is clean. As the ear canal is slightly curved, you have to pull the ear slightly up and backwards before inserting the sensor tip.

The correct position of the thermometer sensor in the ear canal:


- Children **under 1 year of age.**
Lay the child in a horizontal position and turn the head on one side so that the ear is directed upward. Slightly pull the ear back.
 - Children **over 1 year of age.**
Slightly pull the ear up and back to align the ear canal.
4. Carefully insert the sensor tip into the ear canal and press «EAR» button.

If the temperature reaches 37.5°C or above, the backlight turns on and the red LED light is on for 3 seconds.

3. OBJECT TEMPERATURE MEASUREMENT

1. Press «EAR» button to switch on the device.
2. Remove the cap.
3. Simultaneously press and hold for 3 seconds «SET» and «MEM» buttons icon  will appear on the display.
4. Bring the thermometer to the measurement object so that there is a distance of 3-5 cm between the sensor and the object.
5. Press «EAR» button.

When the measurement is completed, you will hear a beep sound, and the measured value will appear on the display.

OPERATION VIA MOBILE APPLICATION

NOTE

Mobile application provides opportunity to:

- save the temperature measurement history;
- save the temperature measurement history for different users;
- take notes on symptoms and general state;
- take notes on medicines taken.

NOTE

The object temperature measurement history is not stored in the application.

Install the application: Android 5.0 and higher, iOS 10.0 and higher.

Option A: Open the PlayMarket or AppStore, enter «AGU» in the search box and download the application.

Option B: For direct link to download the application, scan this QR code.

1. Download and run the application.
2. Register and log in with social media site.
3. Select device to connect.
4. Enter the required user data.
5. Switch on the device and press «Connect» button.

6. Follow the application instructions.

10 DEVICE CLEANING

⚠ ATTENTION

Do not use chemically active detergents.
to clean the device
Disinfecting agents can be used in
the home.

Lens/measuring sensor: In case of contamination of the infrared thermometer sensor lens, the readings begin to change immediately as something interferes with the measurement (usually the readings are too low). To eliminate this problem, it is enough to wipe the sensor with a soft lint-free cloth (or cotton swab) moistened with disinfectant solution. Then it is necessary to wipe the lens with a dry soft lint-free cloth. After the cleaning procedure, the lens surface should be

crystal clear. As the lens surface gets cool when the disinfectant solution evaporates, the thermometer should be left for several minutes at room temperature before use after cleaning.

Thermometer: Clean the thermometer housing with a soft, dry cloth. Never use abrasive cleaners, thinners or benzene. Do not scratch the surface of the sensor lens or display. Avoid dropping your device and do not expose it to extreme temperatures, moisture and direct sunlight.

Store the device in a dry, clean and dust-free place in accordance with the operating conditions.

11 TROUBLESHOOTING

If a problem occurs during operation, first of all, refer to the list of possible failures given below.

LCD display	Cause	Solution
	<ol style="list-style-type: none"> The human body temperature measured is higher than 42°C (107.6°F). The object temperature measured is higher than 100°C (212°F) 	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Service
	<ol style="list-style-type: none"> The human body temperature measured is lower than 35°C (95°F). The object temperature measured is lower than 10°C (50°F) 	
	The operating temperature is out of range between 15°C and 35°C (59°F ~ 95°F)	Operate the thermometer only between the specified temperature ranges

12 STORAGE, TRANSPORTATION AND OPERATION TERMS

- The device should be stored within the temperature range from -25°C to $+55^{\circ}\text{C}$ (-13°F ... $+131^{\circ}\text{F}$) with relative humidity up to 95% (non-condensing).
- The device should be operated within the temperature range from $+15^{\circ}\text{C}$ to $+35^{\circ}\text{C}$ (59°F ~ 95°F) with relative humidity up to 95% (non-condensing).
- The device should be transported within the temperature range from -25°C to $+55^{\circ}\text{C}$ (-13°F ... $+131^{\circ}\text{F}$) with relative humidity up to 85% (non-condensing).
- Do not expose the device to thermal shock.

ATTENTION

After transportation or storage at low temperatures, it is necessary to keep the device at room temperature for at least 2 hours before switching on.

13 DISPOSAL

This symbol on the product or its package means that this product does not fall under the category of domestic waste.

- Proper disposal of the device will prevent adverse environmental and human health effects.
- In order to protect the environment, the device must not be disposed of together with domestic (household) waste. Disposal shall be provided in accordance with local regulations.
- The device must be disposed of in accordance with the EU Directive 2012/19/EU WEEE (Waste Electrical and Electronic Equipment).

If you have any questions, please contact the local public utility responsible for waste disposal.

14 CERTIFICATION

This product conforms to the provisions of the EU Medical Device Directive (93/42/EEC)

EMC Tables

Guidance and manufacturer's declaration – electromagnetic emissions		
Smart infrared thermometer AGU SHE7 is intended for use in the electromagnetic environment specified below. The customer or the user of the AGU SHE7 should assure that it is used in such environment		
Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	Smart infrared thermometer AGU SHE7 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
RF emissions CISPR 11	Class B	Smart infrared thermometer AGU SHE7 is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
Smart infrared thermometer AGU SHE7 is intended for use in the electromagnetic environment specified below. The customer or the user of the AGU SHE7 should assure that it is used in such environment			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact	±6 kV contact	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
	±8 kV air	±8 kV air	

Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
--	-------	-------	--

Radiated RF IEC 61000-4-3	80 MHz to 2.5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of smart infrared thermometer AGU SHE7, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d=1.2\sqrt{P}$ $d=1.2\sqrt{P} \text{ 80 MHz to 800 MHz}$ $d=2.3\sqrt{P} \text{ 800 MHz to 2.5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a) should be less than the compliance level in each frequency range (b) Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
------------------------------	-------------------	-------	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

a. Field strengths from fixed transmitters, such as base stations from radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be adequately predicted by theoretical calculations. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which **AGU SHE7** is used exceeds the applicable RF compliance level above, **AGU SHE7** should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating smart infrared thermometer **AGU SHE7**.

b. Over the frequency range 150kHz to 80MHz, field strengths should be less than [V1] V/m

Recommended separation distances between portable and mobile RF communications equipment and Smart infrared thermometer AGU SHE7

Smart infrared thermometer **AGU SHE7** is intended for use in an electromagnetic environment where radiated RF distances are controlled. The customer or the user of **AGU SHE7** can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and smart infrared thermometer **AGU SHE7** as recommended below according to the maximum output power of the communications equipment

Rated maximum output power of transmitter, W	Separation distance according to frequency of transmitter, m		
	150 kHz – 80 MHz $d=1.2\sqrt{P}$	80 kHz – 800 MHz $d=1.2\sqrt{P}$	800 kHz – 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
10	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people

15 MANUFACTURER'S WARRANTY

The warranty period of this product is 12 months from the date of sale, providing that all operation conditions placed in this instruction manual are strictly observed. The warranty is valid only on presentation of the sales receipt and the warranty card completed by the authorized dealer. The warranty does not apply to consumables, constituent parts and components, such as battery, package, etc.

Warranty and free maintenance service is not provided in case of:

- violation of operating instructions;
- damage caused by intentional or erroneous actions of the consumer due to improper or negligent handling;
- evidence of mechanical impact, dents, cracks, chips, etc. on the thermometer housing, evidence of opening the device, disassembly, unauthorized repair, ingress of moisture, exposure to aggressive agents or any other unauthorized alterations of the device, and in other cases of violation of storage, cleaning, transportation and operation rules specified in the instruction manual.
- ingress of oils, dust, insects, liquids and other foreign objects inside the device.

Follow the instructions precisely to ensure reliable and long-term operation of the device.

In case of abnormal operation of the device, please contact your retailer or Customer Service.

The manufacturer reserves the right to make alterations in the design of the device.

ATTENTION

Do not try to repair the device yourself – this will void the warranty.

The device alterations without the manufacturer's permission will void the warranty.

For repair and maintenance, please contact a specialized after-sales service (see agu-baby.com).

Date of manufacture  and importer information are placed on the individual package.

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AGU[®] is the registered trademark by Montex Swiss AG,
AViTA Corporation, 9F, No. 78, Sec. 1, Kwang-Fu Rd.,
Wellkang Ltd., Suite B, 29 Harley Street, London W1G 9QR, England,